

CardioSecur Pro

22-Lead ECG


Outshining the Gold Standard



Made
in
Germany

360° View of the Heart

CardioSecur Pro is a mobile, 22-lead ECG for medical personnel running on iPhone/iPad. It accomplishes much more than a standard 12-lead ECG. The ECG recordings are comprised of the standard 12-leads plus optionally V7-V9 and VR3-VR9. In addition to complex arrhythmias, CardioSecur makes it possible to recognize infarcts of the anterior, lateral, and posterior walls of the heart. The European Society of Cardiology guidelines recommend examining the posterior and right walls of the heart. CardioSecur is the only system that practically implements this guideline into its measurement process.




The ultra-lightweight cable (50 g) requiring only 4 electrodes and an app is truly one of a kind worldwide. Like the stethoscope, CardioSecur can be used wherever it is needed: **directly with the patient**. Recording a comprehensive ECG has never been so quick and intuitive. The app's communication features are ground-breaking in comparison to those of a conventional ECG. CardioSecur's lead system has been validated in numerous studies. CardioSecur Pro is approved as a class IIa medical product.

Application

The **CardioSecur ECG cable** connects directly to your iPhone/iPad. The optional GDT interface makes it possible to connect to your practice software or Hospital Information System (HIS). CardioSecur also offers the options of ECG interpretation or parameter measurement.

Receive an ECG report in just 3 easy steps:



Thanks to the **manageable size** and easy communication features, you can take CardioSecur with you anytime and anywhere to make a comprehensive diagnosis. With the aid of the zoom function you can analyse the details of individual ECG segments. 10,000 minutes of ECG recordings require 1 GB of storage capacity. These can be archived on your device or database and sent to a colleague as a PDF or printed.

Technical Details

- 12-/22-Lead ECG (I, II, III, aVR, aVL, aVF, V1-V6 plus V7-V9 and VR3-VR9)
- Ultra-lightweight ECG cable (only 50 g) with only 4 electrodes
- ECG recordings of up to 30 minutes each
- Communication via GDT, email, iMessage, iCloud Drive, AirDrop and AirPrint
- ECG System: Bipolar, measurements from 3 channels
- ECG Storage: 10,000 minutes of recording per 1 GB
- Heart rate: 18 to 256 beats per minute
- Bandwidth: 0,05 to 40 Hz, Sampling rate: 250 Hz
- Disposable electrodes: YD43-6-W, Type P-00-S/xx
- CardioSecur is available in the App Store (requires iOS 7.0 or higher)
- Optimised for iPhone, iPad mini, iPad and iPad Pro

CardioSecur uses an effective electrode system based upon the EASI standard^{1,2}. CardioSecur has further developed the EASI standard so that it only requires four electrodes. The positioning of the 4 electrodes capture the entire heart and computes all 22-leads (I, II, III, aVR, aVL, aVF, V1-6, plus V7-V9 and VR3-VR9).

1) Accuracy of the EASI 12-Lead electrocardiogram compared to the standard 12-lead electrocardiogram for diagnosing multiple cardiac abnormalities; J Electrocardiol. 1999; 32 Suppl: 38-47. 2) Comparison of waveforms in conventional 12-lead ECGs and those derived from EASI leads in children; J Electrocardiol. 2003; 36(1): 25-31.

Overview of Benefits

Price Benefits

Instead of a high investment of €3,000 or more for a conventional 12-lead ECG device with set feature usage, CardioSecur offers a flexible feature and price model:

Only 9,95 € / month

+ 199,00 € once
(cancellable monthly)

- Save high one-time costs
- Use a system especially designed for medical personnel
- Option to cancel monthly (without indication of reasons)
- You can also change your individual features (12-leads, 22-leads, ECG interpretation and parameter measurement, etc.) on a monthly basis

Subscription Benefits

- The CardioSecur app and all future updates
- Continuous further development and maintenance of the EU-GDPR compliant medical device platform
- Local archival of your ECG recordings in the app
- Allocate ECGs to each patient and insert comments
- Send or print ECGs directly from your iPhone or iPad (e.g. Email, iChat, Airdrop)
- Option: 22-leads
- Option: ECG interpretation or parameter measurement via the HES algorithm
- Option: data transmission to your Hospital Information System (HIS) via a GDT connection

CardioSecur

Knows Your Heart



CardioSecur is an innovation of Personal MedSystems GmbH
ISO 13485 certified by TÜV SÜD



Personal MedSystems GmbH
Wilhelm-Leuschner-Straße 41
60329 Frankfurt am Main

+49 (0) 69 9072013 0
info@cardiosecur.com

Available at:
www.mobile-ecg.com

Download the CardioSecur Pro app in the App Store:

